

ProStar Plus® XXL Technical Data

Usable Catheter Length	Stent I.D. at Nominal Pressure	Stent Length	Rec. I.D. Introducer Sheath	Rated Burst Pressure
135 cm	[mm]	[mm]	[F]	[bar]
XXL 1000/15-120	10	13	10	8
XXL 1000/20-120	10	18	10	8
XXL 1000/30-120	10	28	10	8
XXL 1000/40-120	10	38	10	8
XXL 1000/50-120	10	48	10	8
XXL 1000/60-120	10	58	10	8
XXL 1200/15-120	12	13	10	8
XXL 1200/20-120	12	18	10	8
XXL 1200/30-120	12	28	10	8
XXL 1200/40-120	12	38	10	8
XXL 1200/50-120	12	48	10	8
XXL 1200/60-120	12	58	10	8
XXL 1400/15-120	14	13	10	8
XXL 1400/20-120	14	18	10	8
XXL 1400/30-120	14	28	10	8
XXL 1400/40-120	14	38	10	8
XXL 1400/50-120	14	48	10	8
XXL 1400/60-120	14	58	10	8
XXL 1600/15-120	16	13	10	8
XXL 1600/20-120	16	18	10	8
XXL 1600/30-120	16	28	10	8
XXL 1600/40-120	16	38	10	8
XXL 1600/50-120	16	48	10	8
XXL 1600/60-120	16	58	10	8
XXL 1800/50-120	18	48	10	8
XXL 1800/60-120	18	58	10	8
XXL 2000/50-120	20	48	10	8
XXL 2000/60-120	20	58	10	8
XXL 2200/50-120	22	48	10	8
XXL 2200/60-120	22	58	10	8
XXL 2400/50-120	24	48	10	8
XXL 2400/60-120	24	58	10	8

ProSTAR PLUS XXL Unmounted Stent:

Order No.	Length of stent
XXL13	13mm
XXL18	18mm
XXL28	28mm
XXL38	38mm
XXL48	48mm
XXL58	58mm

Deployment pressure @ 4 bar



www.vascularconcepts.com

Manufactured By:
 Vascular Concepts Limited
 No. 19, SV Complex
 Bellary Road, Hebbal
 Bangalore 560 024
 INDIA
 Tel: (080) 23438146-49
 Fax: (080) 23439205

