

ProLink SVTM

Balloon Expandable Coronary Small Vessel Stent System

0.070 mm
Articulation Width

0.090 mm
Strut Thickness

0.090 mm
Strut Width

0.30 – 0.35 bar
Radial Strength

VASCULAR
c o n c e p t s

Manufactured By:
Vascular Concepts Ltd.
No. 19, SV Complex
Bellary Road, Hebbal
Bangalore 560 024 India
Tel: (080) 23438146 / 47 / 48 / 49
Fax: (080) 23439205

ProLink SV Coronary Stent System Technical Data

BALLOON [mm]	2.00	2.25	2.50
PRESSURE [Mpa/bar]			
0.8/8	2.00	2.25	2.50
0.9/9	2.04	2.28	2.53
1.0/10	2.08	2.30	2.56
1.1/11	2.10	2.33	2.59
1.2/12	2.13	2.36	2.62
1.3/13	2.15	2.39	2.65
1.4/14	2.19	2.42	2.68
1.5/15	2.21	2.46	2.71
1.6/16	2.23	2.50	2.74
1.7/17	2.25	2.54	2.78
1.8/18	2.28	2.60	2.82
1.9/19	2.31	2.66	

NOMINAL PRESSURE AT 8 BAR (0.8MPa)

RATED BURST PRESSURE

ProLink SV Coronary Stent System Ordering Information

LENGTH	BALLOON DIAMETER		
	2.00 mm	2.25 mm	2.50 mm
10 mm	PLSV2010	PLSV22510	PLSV2510
15 mm	PLSV2015	PLSV22515	PLSV2515
20 mm	PLSV2020	PLSV22520	PLSV2520
25 mm	PLSV2025	PLSV22525	PLSV2525
30 mm	PLSV2030	PLSV22530	PLSV2530