


Speed+: Fast Exchange Semi-Compliant Angioplasty Balloon

Redefining Engineering for Outstanding Trackability



Speed+ Technical Data

PRESSURE [MPa/atm] BALLOON [mm]	DIAMETER (mm)							
	1.5	2.00	2.25	2.50	3.00	3.25	3.50	4.00
0.8/8	1.50	2.00	2.25	2.50	3.00	3.25	3.50	4.00
0.9/9	1.53	2.04	2.28	2.53	3.03	3.28	3.55	4.04
1.0/10	1.56	2.08	2.30	2.56	3.06	3.32	3.59	4.09
1.1/11	1.58	2.10	2.33	2.59	3.09	3.35	3.62	4.12
1.2/12	1.61	2.13	2.36	2.62	3.12	3.38	3.65	4.15
1.3/13	1.64	2.15	2.39	2.65	3.14	3.41	3.69	4.19
1.4/14	1.67	2.19	2.42	2.68	3.17	3.44	3.72	4.23
1.5/15	1.70	2.21	2.46	2.71	3.20	3.47	3.76	4.27
1.6/16	1.73	2.23	2.50	2.74	3.22	3.50	3.80	4.31
1.7/17	1.78	2.25	2.54	2.78	3.25	3.53	3.84	4.36
1.8/18	1.81	2.28	2.60	2.82	3.28	3.56	3.89	4.40
1.9/19	1.85	2.31	2.66	2.86	3.32	3.59	3.94	4.45
2.0/20	1.93	2.34	2.75	2.92	3.35	3.63	4.00	4.50

Nominal Pressure at 8 atm [0.8 MPa] RATED BURST PRESSURE

Speed+ Ordering Information

BALLOON LENGTH [mm]	DIAMETER (mm)							
	1.50	2.00	2.25	2.50	3.00	3.25	3.50	4.00
8	SP1508	SP2008	SP22508	SP2508	SP3008	SP32508	SP3508	SP4008
10	SP1510	SP2010	SP22510	SP2510	SP3010	SP32510	SP3510	SP4010
12	SP1512	SP2012	SP22512	SP2512	SP3012	SP32512	SP3512	SP4012
15	SP1515	SP2015	SP22515	SP2515	SP3015	SP32515	SP3515	SP4015
20	SP1520	SP2020	SP22520	SP2520	SP3020	SP32520	SP3520	SP4020
25	SP1525	SP2025	SP22525	SP2525	SP3025	SP32525	SP3525	SP4025
30	SP1530	SP2030	SP22530	SP2530	SP3030	SP32530	SP3530	SP4030



www.vascularconcepts.com

Vascular Concepts Limited
No. 19, SV Complex
Bellary Road, Hebbal
Bangalore 560 024
INDIA
Tel: (080) 23438146-49
Fax: (080) 23439205

