


ProTAXX

Paclitaxel Eluting Coronary Stent System

**VASCULAR**
c o n c e p t s
...The Innovations Company

ProTAXX: Paclitaxel Eluting Coronary Stent System

Finally a drug-eluting stent that leaves nothing behind.



ProTAXX Compliance Data

STENT I.D (mm) Pressure Mpa/atm	Compliance (mm)									
	2.50	2.75	3.00	3.25	3.50	4.00	4.50	5.00	5.50	6.00
0.6/6	2.42	2.68	2.92	3.18	3.40	3.92	4.38	4.90	5.4	5.88
0.7/7	2.46	2.71	2.96	3.22	3.45	3.96	4.44	4.95	5.45	5.94
0.8/8	2.5	2.75	3.00	3.25	3.50	4.00	4.50	5.00	5.50	6.00
0.9/9	2.53	2.78	3.03	3.28	3.55	4.04	4.56	5.05	5.55	6.06
1.0/10	2.56	2.81	3.06	3.32	3.59	4.09	4.60	5.09	5.60	6.11
1.1/11	2.59	2.84	3.09	3.35	3.62	4.12	4.64	5.13	5.65	6.16
1.2/12	2.62	2.86	3.12	3.38	3.65	4.15	4.68	5.16	5.69	6.22
1.3/13	2.65	2.89	3.14	3.41	3.69	4.19	4.72	5.20	5.73	6.27
1.4/14	2.68	2.91	3.17	3.44	3.72	4.23	4.76	5.23	5.77	6.33
1.5/15	2.71	2.94	3.20	3.47	3.76	4.27	4.81	5.27	5.82	6.39
1.6/16	2.74	2.97	3.22	3.50	3.80	4.31	4.86	5.31	5.87	6.45
1.7/17	2.78	3.00	3.25	3.53	3.84	4.36	4.92	5.35	5.93	6.52
1.8/18	2.82	3.03	3.28	3.56	3.89	4.40	4.98	5.40	5.99	6.60
1.9/19	2.86	3.07	3.32	3.59	3.94	4.45	5.04	5.45	6.05	6.68
2.0/20	2.92	3.11	3.35	3.63	4.00	4.50	5.12	5.50	6.11	6.75

Nominal Pressure at 8 atm [0.8 MPa]

RATED BURST PRESSURE

ProTAXX Ordering Information

STENT I.D.	STENT LENGTH							
	8mm	10mm	13mm	16mm	18mm	23mm	28mm	33mm
2.50mm	PT2508	PT2510	PT2513	PT2516	PT2518	PT2523	PT2528	PT2533
2.75mm	PT27508	PT27510	PT27513	PT27516	PT27518	PT27523	PT27528	PT32533
3.00mm	PT3008	PT3010	PT3010	PT3016	PT3018	PT3023	PT3028	PT3033
3.25mm	PT32508	PT32510	PT32513	PT32516	PT32518	PT32523	PT32528	PT32533
3.50mm	PT3508	PT3510	PT3513	PT3516	PT3518	PT3523	PT3528	PT3533
4.00mm	PT4008	PT4010	PT4013	PT4016	PT4018	PT4023	PT4028	PT4033
4.50mm	—	—	PT4513	PT4516	PT4518	PT4523	PT4528	PT4533
5.00mm	—	—	PT5013	PT5016	PT5018	PT5023	PT5028	PT5033
5.50mm	—	—	PT5513	PT5516	PT5518	PT5523	PT5528	PT5533
6.00mm	—	—	PT6013	PT6016	PT6018	PT6023	PT6028	PT6033



www.vascularconcepts.com

Manufactured By:
Vascular Concepts Limited
2 Deans Hall Business Park
Oak Road, Little Maplestead
Halstead, Essex C09 2RT
United Kingdom
Tel: 01787 477199
Fax: 01787 477510



Vascular Concepts Limited
No. 19, SV Complex (1st Floor)
Bellary Road, Hebbal
Bangalore 560 024
INDIA
Tel: (080) 23438146-49
Fax: (080) 23439205